

PHYS 461 — Modern Physics Lab

Cavendish Experiment

Micah Hillman

Cordney Nash

Report December 6, 2024

Performed November 2024

Introduction

The Cavendish experiment, the namesake experiment of 18th century physicist Henry Cavendish, successfully measured the density of Earth, unintentionally allowing for the first measurement of Newton’s universal gravitational constant G — calculated by later physicists [2]. The experiment uses a hermetically sealed torsion balance system to track the slow oscillation induced by bringing large masses (in our case, lead) in close proximity to a set of test masses connected to a torsional balance; by these means, the masses are made to interact gravitationally while minimizing the influence of extraneous forces [1]. Bringing the test masses in close proximity changes the torsional pendulum’s equilibrium position, to which the ensuing oscillations gradually decay (with a regular period — T — determined by the balance’s moment of inertia and torsional constant) [1]. By recording these oscillations and measuring the gravitationally-modulated torsional balance’s changing equilibrium position, the universal gravitational constant G may be calculated [1].

Theory

At its core, the Cavendish experiment relies on forces described by Newtonian universal gravitation, famously expressed in the equation

$$F_g = G \frac{mM}{r^2}.$$

When large lead masses (M) are positioned in close proximity to the torsional balance’s smaller masses (m), the equilibrium position of the torsional balance is suddenly altered, leading to slow *oscillatory behavior* about a new equilibrium — the period (T) of which depends on the balance’s torsional constant (κ) and moment of inertia ($I = 2md^2$):

$$T = 2\pi \sqrt{\frac{2md^2}{\kappa}}. \tag{1}$$

The torsional constant (κ) can be expressed as

$$\kappa = \frac{4GmMd}{r^2 \Delta\theta}, \tag{2}$$


so we can combine these two expressions (1, 2) and isolate G in terms of observables and experimental parameters [1]:

$$G = \frac{2\pi^2 r^2 d\Delta\theta}{MT^2}. \quad (3)$$


Methods

22kg lead sphere

Results



(a). Birds-eye illustration of the sealed torsional balance system and associated weights. Important components include both the small and large lead masses m and M , a reflective mirror attached to the central portion of the torsional balance beam, and a laser light which is reflected across the room to a measuring stick.



(b). Illustration depicting measurement. A camera (center) and meter stick (left, at a known distance) were used to capture the position of laser light reflected from the mirror affixed to central axis of the torsional balance (right, and shown in figure 1a).

Figure 1. Two sketched views of the Cavendish apparatus. Figure 1a shows components associated with torsional balance system, while figure 1b depicts components associated with measurement.

Bibliography

- [1] BlueMarbleScience. *Calculating the Universal Gravitational Constant*. <https://www.youtube.com/watch?v=W2uPNg1BUSI>. Retrieved December 4, 2024. Dec. 2020.
- [2] American Physical Society. *1798: Cavendish Weighs the World*. <https://www.aps.org/apsnews/2008/06/1798-cavendish-weighs-world>. Retrieved December 4, 2024. June 2008.